

PATIENT.

COMMENCING DIALYSIS VIA CYCLING MACHINE

- REQUIREMENTS: - Cycling machine  
pack of lines  
drainage line or bag  
dialysis fluid  
~~green shield~~ *connection shield*  
trolley wipes  
Heparin,  
2 ml syringe  
needles x 2  
sterets x 2  
packet gauze  
exit site dressing  
Betadine scrub

1. Wash hands and collect equipment.
2. Clean machine.
3. Hang up dialysis fluid onto machine.
4. Line machine as directed on card and confirm programme.
5. Short hand wash and clean trolley.
6. Set out requirements for exit care and addition of Heparin. *Connection shield*
7. Expose dialysis catheter. Remove exit dressing, inspect, clean and redress site. *Spring Green Blue redk Heparin scrub*
8. Wash hands thoroughly for 1 minute, rinse and repeat wash.
9. Draw up Heparin and add to "heater" bag.
10. Connect lines to bags - connecting the "patient line" to heater bag first. (If additional bags required during treatment, then place green shield onto 1 bag at this stage).
11. Prime patient line. *Break Green seal on heater bag & open while clamp on patient line*
12. Hibisol hands and open ~~shield~~ *connection* pack. *x 2*
13. Apply ~~shield~~ *connection* to patient line.
14. Remove disconnect cap or bag from patient and attach to patient line at machine.
15. Press "start/continue"..
16. Open all clamps (remember roller clamp at patient). Press "start/~~continue~~" and commence treatment.

Purging Tubing  
Drawn 0 of 08